

Rat anti-Mouse IgM, clone LO-MM-9, Biotin (Monoclonal)

Clone no. LO-MM-9

MONOSAN

Product name	Rat anti-Mouse IgM, clone LO-MM-9, Biotin (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1k
Clonality	Monoclonal
Clone number	LO-MM-9
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgM, clone LO-MM-9, Biotin (Monoclonal)

Clone no. LO-MM-9

MONOSAN

Additional info

Mouse heavy chain of immunoglobulin (determined by immunodot). Avidity on IgM: 3.6×10^9 M⁻¹. No crossreactivity with other animals. Available as unconjugated MON5057, PE conjugated MON 5057P, Biotin conjugated MON 5057B, FITC conjugated MON5057F

References

1. Bazin et al. Pergamon Press Oxford and NY 1982; 29:615-618
2. Cormont et al. Meth in Enzymol. 1986;121:622-631
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES