

Rat anti-Mouse IgM, clone LO-MM-8, HRP (Monoclonal)

Clone no. LO-MM-8

MONOSAN

Product name	Rat anti-Mouse IgM, clone LO-MM-8, HRP (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2ak
Clonality	Monoclonal
Clone number	LO-MM-8
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgM, clone LO-MM-8, HRP (Monoclonal)

Clone no. LO-MM-8

MONOSAN

Additional info

Mouse heavy chain of immunoglobulin (determined by immunodot). Avidity on IgM: 7×10^8 M⁻¹. No crossreactivity with other animals. Available in unconjugated MON5057, HRP conjugated MON5057P

References

1. Bazin et al. Prot Biol Fluis 1982; 615-618
2. Cormont et al. Meth in Enzymol. 1986;121:622-631
3. Savignac et al. Journal of Immunol. 2010
4. Huang et al. Journal of Immunol. 2011
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES