

Rat anti-Mouse IgG2a, clone LO-MG-2a9, HRP (Monoclonal)

Clone no. LO-MG-2a9

MONOSAN

Product name	Rat anti-Mouse IgG2a, clone LO-MG-2a9, HRP (Monoclonal)
Host	Rat
Applications	IHC-fr, ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1k
Clonality	Monoclonal
Clone number	LO-MG-2a9
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS pH 7.2 / 50% (v/v) Glycerol
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse IgG2a, clone LO-MG-2a9, HRP (Monoclonal)

Clone no. LO-MG-2a9

MONOSAN

Additional info

Mouse gamma 2a heavy chain of immunoglobulin (determined by immunodot). No crossreactivity with other animals. Useful in immunoassays, such as ELISA's (good binding on plastic), immunohistochemistry on frozen sections and for affinity chromatography of murine IgG1. Available unconjugated MON5068, HRP conjugated 5068H

References

1. Bazin, et al., J.Immunology Meth. 1984;71, 9-16
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES