

Rat anti-Mouse Interleukin-1 receptor/P80 clone Reg21 (Monoclonal)

Clone no. Reg21

MONOSAN

---

Product name	Rat anti-Mouse Interleukin-1 receptor/P80 clone Reg21 (Monoclonal)
Host	Rat
Applications	FUNC,ELISA
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	Reg21
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse Interleukin-1 receptor/P80 clone Reg21 (Monoclonal)

Clone no. Reg21

MONOSAN

---

**Additional info**

The monoclonal antibody binds to soluble mouse interleukin-1 receptor type I. The IL-1 system includes two agonists (IL-1alpha and IL-1beta), converting enzymes, antagonists, two receptors (IL-1RI and IL-1RII) and the IL-1 receptor accessory protein. Interleukin-1 signal is transduced through the type I receptor.

**References**

1. Kojouharoff, G et al. Clin Exp Immunol 1997; 107: 353
2. Stern, A et al Arch Biochem Biophys 1989, 274: 26
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES