

Rat anti Mouse-IgG2a Kappa, clone LO-MK-1, FITC (Monoclonal)

Clone no. LO-MK-1

MONOSAN

Product name	Rat anti Mouse-IgG2a Kappa, clone LO-MK-1, FITC (Monoclonal)
Host	Rat
Applications	IHC, ELISA, FC, IF, WB
Species reactivity	mouse
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2ak
Clonality	Monoclonal
Clone number	LO-MK-1
Size	1 ml
Concentration	1 mg/ml
Format	-
Storage buffer	PBS with 0.1% sodium azide
Storage until expiry date	-20°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti Mouse-IgG2a Kappa, clone LO-MK-1, FITC (Monoclonal)

Clone no. LO-MK-1

MONOSAN

Additional info

Mouse gamma light chain of Immunoglobulin (determined by immunodot) and exhibits no reactivity to mouse lambda light chain. No detectable reactivity with human IgG as tested by ELISA. Available in unconjugated MON5083, FITC conjugated MON5083F, HRP conjugated MON5063H and Biotin conjugated MON5063B.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES