Product datasheet

MON5164-1



MONOSAN

Mouse anti Pan Cytokeratin, clone BS5 (Monoclonal)Clone no.BS5

Product name	Mouse anti Pan Cytokeratin, clone BS5 (Monoclonal)	
Host	Mouse	
Applications	IHC-P (1:100-1:400), IHC-Fr	
Species reactivity	human, sheep, pig, dog, mouse	
Conjugate	-	
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers	
lsotype	lgG1	
Clonality	Monoclonal	
Clone number	BS5	
Size	1 ml	
Concentration	n/a	
Format	-	
Storage buffer	TRIS with 0.03% sodium azide, pH7,2	
Storage until expiry date	2-8°C	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON5164-1



Mouse anti Pan Cytokeratin, clone BS5 (Monoclonal)Clone no.BS5

MONOSAN

Additional info

Cytokeratins are classified into one of two classes, type I (acidic polypeptides) and type II (basic polypeptides). Cytokeratins play a critical role in differentiation and tissue specialization and function to maintain the overall structural integrity of epithelial cells. Cytokeratins have been found to be useful markers of tissue differentiation which is directly applicable to the characterization of malignant carcinomas.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com