Product datasheet MON6030



Mouse anti-Double stranded DNA, clone AE-2 (Monoclonal)

Clone no. AE-2 MONOSAN

Product nameMouse anti-Double stranded DNA, clone AE-2 (Monoclonal)

Host Mouse

Applications FC, IHC-fr, IHC-P, IF

Species reactivity human, mouse, rat, rabbit, all eukaryotes

Conjugate -

Immunogen nuclei of Raji Burkitt's lymphoma cells

lsotype lgG3

Clonality Monoclonal

Clone number AE-2

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON6030



Mouse anti-Double stranded DNA, clone AE-2 (Monoclonal)

Clone no. AE-2 MONOSAN

Additional info

AE-2 recognizes double stranded DNA expressed in the nucleus of

eukaryotes

References 1. Epstein, AL et al. In Progress on nonhistone protein research. 1: 117-137 (198)

2 -

3.

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES