

Mouse anti-Matrix Metalloproteinase-1, clone 3B6 (Monoclonal)

Clone no. 3B6

MONOSAN

Product name	Mouse anti-Matrix Metalloproteinase-1, clone 3B6 (Monoclonal)
Host	Mouse
Applications	IHC-P, WB
Species reactivity	Human
Conjugate	-
Immunogen	human MMP-1 peptide: CSSFGFPRTVKH
Isotype	IgG1
Clonality	Monoclonal
Clone number	3B6
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with 0.05% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Matrix Metalloproteinase-1, clone 3B6 (Monoclonal)

Clone no. 3B6

MONOSAN

Additional info

This monoclonal antibody binds to human MMP-1. Reactivity in other species has not been determined. The antibody was tested for cross-reactivity with MMP-2, MMP-3, and MMP-9 and did not cross-react.

References

1. Murray GI et al. J Pathol 1998;185:256-261
2. Murray GI et al. GUT 1998;43:791-797
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES