Product datasheet

MON6058



MONOSAN

Mouse anti-MUC5AC, clone 45M1 (Monoclonal)Clone no.45M1

Product name	Mouse anti-MUC5AC, clone 45M1 (Monoclonal)	
Host	Mouse	
Applications	ELISA, IHC-fr, IHC-P, WB	
Species reactivity	human, mouse, rat, rabbit, monkey, chicken, cow, cat, pig, hedgehog	
Conjugate	-	
Immunogen	mucin isolated from an ovarian cyst fluid	
lsotype	lgG1	
Clonality	Monoclonal	
Clone number	45M1	
Size	100 ug	
Concentration	100 ug/ml	
Format	-	
Storage buffer	PBS with 0.02% sodium azide	
Storage until expiry date	2-8°C	

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON6058

MONOSAN

Mouse anti-MUC5AC, clone 45M1 (Monoclonal) Clone no. 45M1

MONOSAN

Additional info

45M1 recognizes the peptide core of gastric mucin M1 (now: MUC5AC), and more specifically with the 'h' epitope amongst the a, b, c, d, e, f, g and h protein core epitopes defined by Bara for M1. 45M1 and 2-12M1 both specifically react with epitopes located in the C-terminal cysteine rich part of the peptide core of MUC5AC. MUC5AC is present in primary ovarian mucinous cancer and gastric cancer, but usually absent in colorectal adenocarcinoma, thus showing an expression pattern opposite to MUC2. Anti-MUC5AC may be useful for differential identification of primary mucinous ovarian tumors from colon adenocarcinoma metastatic to the ovary. MUC5AC antibodies may also be useful for identification pancreatic carcinoma and pre-cancerous changes vs. normal pancreas.

References	1.	Bara, J. et al., Int. J. Cancer 47: 304-310 (1991)
	2	Bara, J. et al., J. Immunol. Methods 149: 105-113 (1992)

- 3. -4. -
- 4. -5. -
- 5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com