

Mouse anti-MMP-3, clone 1B4 (Monoclonal)

Clone no. 1B4

MONOSAN

Product name	Mouse anti-MMP-3, clone 1B4 (Monoclonal)
Host	Mouse
Applications	IHC-P, IHC-fr, WB
Species reactivity	Human
Conjugate	-
Immunogen	human MMP-3 peptide: KSLRKLEPELH
Isotype	IgG1
Clonality	Monoclonal
Clone number	1B4
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with 0.05% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-MMP-3, clone 1B4 (Monoclonal)

Clone no. 1B4

MONOSAN

Additional info

The antibody binds to human MMP-3. Reactivity in other species has not been determined. The antibody was tested for cross-reactivity with MMP-1, MMP-2 and MMP-9 and did not cross-react.

References

1. Murray GI et al. GUT 1998;43:791-797
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES