

## Mouse anti-TIMP-1, clone 2A5 (Monoclonal)

Clone no. 2A5

MONOSAN

---

Product name	Mouse anti-TIMP-1, clone 2A5 (Monoclonal)
Host	Mouse
Applications	IHC-P, WB
Species reactivity	Human
Conjugate	-
Immunogen	human TIMP-1 peptide: FQALGDAADIR
Isotype	IgG1
Clonality	Monoclonal
Clone number	2A5
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with 0.05% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Mouse anti-TIMP-1, clone 2A5 (Monoclonal)

Clone no. 2A5

MONOSAN

**Additional info**

This monoclonal antibody binds to human TIMP-1. Reactivity in other species has not been determined. The 11-mer peptide used for the immunization differs on at least 3 amino acids with the sequence of TIMP-1 of rat, mouse, and rabbit. Therefore it is not expected to be effective in these species as well. The antibody was tested for cross-reactivity with TIMP-2 and did not cross-react.

**References**

1. Murray GI et al. GUT 1998;43:791-797
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES