

## Mouse anti-TIMP-2, clone 3A4 (Monoclonal)

Clone no. 3A4

MONOSAN

---

Product name	Mouse anti-TIMP-2, clone 3A4 (Monoclonal)
Host	Mouse
Applications	IHC-P, IHC-fr, WB
Species reactivity	Human
Conjugate	-
Immunogen	human TIMP-2 peptide: DSGNDIYGNPI
Isotype	IgG1
Clonality	Monoclonal
Clone number	3A4
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with 0.05% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Mouse anti-TIMP-2, clone 3A4 (Monoclonal)

Clone no. 3A4

MONOSAN

**Additional info**

The antibody binds to human TIMP-2. Reactivity in other species has not been determined, but the 11-mer peptide used for the immunization shows a 100% match with rat, mouse, and rabbit TIMP-2. The antibody was tested for cross-reactivity with TIMP-1 and did not cross-react.

**References**

1. Murray GI et al. GUT 1998;43:791-797
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES