

Mouse anti-Rat CD8, clone OX8-2a (Monoclonal)

Clone no. OX8-2a

MONOSAN

---

Product name	Mouse anti-Rat CD8, clone OX8-2a (Monoclonal)
Host	Mouse
Applications	IHC-P, IHC-fr
Species reactivity	Rat
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2a
Clonality	Monoclonal
Clone number	OX8-2a
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with 0.05% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Rat CD8, clone OX8-2a (Monoclonal)

Clone no. OX8-2a

MONOSAN

**Additional info**

The antibody binds to a rat strain-independent epitope of the  $\alpha$ -chain of CD8, present on the CD8-positive subsets of T cells and NK cells.

**References**

1. Kuppen PJ et al. Cancer Immunol 1993;36:403-408
2. Brideau RJ et al. EurJImmunol 1980;10:609-615
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES