

Mouse anti-Rat CD45, clone ANK74 (Monoclonal)

Clone no. ANK74

MONOSAN

Product name	Mouse anti-Rat CD45, clone ANK74 (Monoclonal)
Host	Mouse
Applications	IHC-fr, FC
Species reactivity	Rat
Conjugate	-
Immunogen	IL-2-activated cultured NK cells of Wag rats
Isotype	IgG1
Clonality	Monoclonal
Clone number	ANK74
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with 0.05% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Rat CD45, clone ANK74 (Monoclonal)

Clone no. ANK74

MONOSAN

Additional info

This is a rat strain-independent antibody to CD45. This monoclonal antibody recognizes a common epitope of rat CD45, present on all rat leucocytes. The ANK74 monoclonal antibody was generated by immunizing mice with IL-2-activated cultured NK cells of Wag rats.

References

1. Giezeman-Smits KM et al. J Immunol 1999;163:71-76
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES