

## Mouse anti-Hepatitis B virus, clone HB3 (Monoclonal)

Clone no. HB3

MONOSAN

---

Product name	Mouse anti-Hepatitis B virus, clone HB3 (Monoclonal)
Host	Mouse
Applications	ICC, ELISA
Species reactivity	Hepatitis B Virus
Conjugate	-
Immunogen	Purified HBsAg from human plasma.
Isotype	IgG1
Clonality	Monoclonal
Clone number	HB3
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Hepatitis B virus, clone HB3 (Monoclonal)

Clone no. HB3

MONOSAN

---

**Additional info**

HBsAg (hepatitis B surface antigen) is a marker of infectivity. Its presence indicates either acute or chronic HBV infection.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES