

Rat anti-Laminin, clone MEC5 (Monoclonal)

Clone no. MEC5

MONOSAN

---

Product name	Rat anti-Laminin, clone MEC5 (Monoclonal)
Host	Rat
Applications	IHC-fr
Species reactivity	Rat, Human, Mouse
Conjugate	-
Immunogen	auto-antibody
Isotype	IgG1
Clonality	Monoclonal
Clone number	MEC5
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Culture supernatant with 0.05% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Laminin, clone MEC5 (Monoclonal)

Clone no. MEC5

MONOSAN

---

**Additional info**

The antibody is directed against the P1 fragment of laminin. They are species-independent as it has been shown that they bind to laminin of rat, mouse and humans. MEC5 is an auto-antibody derived from DZB rats that had developed membranous glomerulopathy upon exposure to mercuric chloride

**References**

1. Aten J, et al. Am J Pathol; 146: 1467-1480
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES