

Mouse anti-Influenza A Nucleoprotein, clone AA5H (Monoclonal)

Clone no. AA5H

MONOSAN

Product name	Mouse anti-Influenza A Nucleoprotein, clone AA5H (Monoclonal)
Host	Mouse
Applications	WB,IHC-P,IF
Species reactivity	viral
Conjugate	Purified
Immunogen	Influenza A / Puerto Rico / 8 / 34 (H1N1) and A/Bangkok / 1 / 79 (H3N2) viruses.
Isotype	IgG2a
Clonality	Monoclonal
Clone number	AA5H
Size	1 mg
Concentration	1.0 mg/ml
Format	-
Storage buffer	PBS
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Influenza A Nucleoprotein, clone AA5H (Monoclonal)

Clone no. AA5H

MONOSAN

Additional info

Mouse anti Influenza A Nucleoprotein antibody, clone AA5H recognizes an epitope within Influenza virus A nucleoprotein. Mouse anti Influenza A Nucleoprotein antibody, clone AA5H can be used in influenza A IFA typing in conjunction with MCA401 (clone GA2B).

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES