

Mouse anti-CA15.3, clone DBMM1 (Monoclonal)

Clone no. DBMM1

MONOSAN

---

Product name	Mouse anti-CA15.3, clone DBMM1 (Monoclonal)
Host	Mouse
Applications	IHC-P
Species reactivity	human
Conjugate	-
Immunogen	Human Breast cancer cell line ZR-75-1
Isotype	IgG1k
Clonality	Monoclonal
Clone number	DBMM1
Size	1 ml
Concentration	lot specific
Format	-
Storage buffer	Tissue culture supernatant with sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CA15.3, clone DBMM1 (Monoclonal)

Clone no. DBMM1

MONOSAN

**Additional info**

This antibody is specific to epithelial membrane antigen (EMA), CA 15-3, or polymorphic epithelial mucin. This antibody stains an underglycosylated MUC1 often present on carcinoma cells.

**References**

1. Hilkens et al. 1984;34:197-206
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES