

Mouse anti-Macrophage Marker, clone 3A5 (Monoclonal)

Clone no. 3A5

MONOSAN

Product name	Mouse anti-Macrophage Marker, clone 3A5 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P, IF (1:20-1:50)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG2b
Clonality	Monoclonal
Clone number	3A5
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	ascites with 1% BSA and 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Macrophage Marker, clone 3A5 (Monoclonal)

Clone no. 3A5

MONOSAN

Additional info

This antibody stains macrophages and histocytic cells in a wide variety of normal and diseased tissue. Dendritic cells, myeloid cells and glia cells do not react with this monoclonal antibody. Monocytes stain weakly with 3A5, also after external activation.

References

1. Jaspers et al. J Clin Pathol 1994; 47: 248-252
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES