

Mouse anti-Macrophages, clone LN-5 (Monoclonal)

Clone no. LN-5

MONOSAN

Product name	Mouse anti-Macrophages, clone LN-5 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IHC-P, IF
Species reactivity	human
Conjugate	-
Immunogen	deparaffinized cell extracts of peripheral blood mononuclear cells
Isotype	IgM
Clonality	Monoclonal
Clone number	LN-5
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Macrophages, clone LN-5 (Monoclonal)

Clone no. LN-5

MONOSAN

Additional info

LN-5 reacts with human macrophages and displays lymph node germinal center and mantle zone B cell reactivity. It reacts with interdigitating reticulum cells, with tingible body and sinus histiocytes. It further reacts with certain tumor cells and also with normal nonlymphoid tissue like chief cells of the stomach and spermatogonia. LN-5 is negative on Hodgkin's disease and non-Hodgkin's lymphomas.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES