

Mouse anti-Macrophage, clone EP-3 (Monoclonal)

Clone no. EP-3

MONOSAN

Product name	Mouse anti-Macrophage, clone EP-3 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P, IF
Species reactivity	human
Conjugate	-
Immunogen	human thymus
Isotype	IgM
Clonality	Monoclonal
Clone number	EP-3
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Macrophage, clone EP-3 (Monoclonal)

Clone no. EP-3

MONOSAN

Additional info

Monoclonal antibody EP-3 recognizes an antigen associated with the cytoplasm of human macrophages resident in the thymus and other lymphoid tissues. Monoclonal antibody EP-3 produces a strong cytoplasmic staining pattern of cortical thymic macrophages in formalin fixed, paraffin embedded tissues specimens. It may therefore be used as a marker of this cell type and of tumors derived from the macrophage.

References

1. Fleming MG, et al, J Cutan Pathol. 17(2): 77-81 (1990)
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES