Product datasheet MON7055



## Rat anti-Mouse LBP, clone M330-19 (Monoclonal)

Clone no. M330-19 MONOSAN

Product name Rat anti-Mouse LBP, clone M330-19 (Monoclonal)

**Host** Rat

**Applications** FUNC,ELISA

Species reactivity mouse

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

**Isotype** IgG2a

**Clonality** Monoclonal

Clone number M330-19

Size 1 ml

Concentration 100 ug/ ml

Format -

Storage buffer PBS with 0.1% BSA and 0.02% sodium azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON7055



Rat anti-Mouse LBP, clone M330-19 (Monoclonal)

Clone no. M330-19 MONOSAN

## Additional info

The monoclonal antibody M330-19 reacts highly specific with mouse natural and recombinant LBP. The antibody is a type I antibody blocking the LPS binding to LBP. LPS binding protein (LBP) is an approximately 60 kDa acute phase protein that is produced by hepatocytes. This protein strongly binds to LPS and has been shown to play an important role in the handling of LPS by the host. A number of functions of LBP have been reported. First, LBP transfers LPS to the LPS receptor CD14 on mononuclear phagocytes, leading to an 100-1,000-fold increased sensitivity of the cells to LPS. Furthermore, LBP can enhance the response of CD14 negative cells by acceleration of LPS binding to soluble CD14, a complex that stimulates these cells. Next, LBP transfers LPS into High Density Lipoprotein (HDL), which effectively neutralizes its biological potency. LBP was demonstrated to protect mice from septic shock caused by LPS or gram negative bacteria.

References 1. Le Roy; D et al. J Immunol 1999; 162: 7454

2 Le Roy, D et al J Immunol 2001, 169: 2759

3. Heumann; D et al. Inf and Immun 2001; 69:378

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES