### Product datasheet

MON7088



MONOSAN

Mouse anti Pepsinogen I, clone 4C6.1 (Monoclonal) Clone no. 4C6.1

Product name	Mouse anti Pepsinogen I, clone 4C6.1 (Monoclonal)	
Host	Mouse	
Applications	IHC-fr, IHC-P, WB, ELISA	
Species reactivity	human	
Conjugate	-	
Immunogen	Pepsinogen I purified from human tissue	
lsotype	lgG1	
Clonality	Monoclonal	
Clone number	4C6.1	
Size	1 ml	
Concentration	100 ug/ ml	
Format	-	
Storage buffer	PBS with 1% BSA & 0.09% sodium azide	
Storage until expiry date	2-8°C	

# FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

### Product datasheet

MON7088

# MONOSAN

Mouse anti Pepsinogen I, clone 4C6.1 (Monoclonal) Clone no. 4C6.1

MONOSAN

### Additional info

The antibody is specific to human pepsinogen I and has no crossreactivity to human pepsinogen II. Pepsinogen I is a group of precursor molecules for pepsin. These proteins are solely synthetized and secreted into gastric lumen by chief (pepsin) cells and mucous neck cells in the gastric corpus (oxyntic mucosa). In atrophic corpus gastritis these cells disappear resulting in a decrease of the serum level of pepsinogen I and in a reduction of the number of pepsinogen I positive cells in gastric biopsies. The presence of positive immunostaining for pepsinogen I is a highly reliable sign for the acidsecreting oxyntic glands. In gastric heterotopia of the duodenal bulb, but not in gastric metaplasia, the oxyntic-type glands give a positive immunohistochemical reaction for pepsinogen I.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com