

## Mouse anti-Cellular Fibronectin, clone DH1 (Monoclonal)

Clone no. DH1

MONOSAN

---

Product name	Mouse anti-Cellular Fibronectin, clone DH1 (Monoclonal)
Host	Mouse
Applications	IHC-fr, WB, ELISA
Species reactivity	human, chicken, rat, rabbit, guinea pig
Conjugate	-
Immunogen	Purified fibronectin from A8387 fibrosarcoma cells
Isotype	IgG1
Clonality	Monoclonal
Clone number	DH1
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 1% BSA & 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Cellular Fibronectin, clone DH1 (Monoclonal)

Clone no. DH1

MONOSAN

---

**Additional info**

The antibody is specific to extradomain A (EDA) sequence of a cellular fibronectin and recognizes thus only the cellular fibronectin. It has been shown that it specifically block chondrocyte condensation in chicken embryos

**References**

1. Vartio et al. J Cell Sci 1987;88:419-430
2. Virtanen et al. Histochem 1988;90:25-30
3. Gould et al. Cell Diff Devel 1990;32:409-416
4. Laitinen et al. Lab Invest 1991;64:492-498
5. Ylatupa et al. Clin Chim Acta 1995;234:79-90

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES