

Mouse anti-Beta1-Integrin, clone DF5 (Monoclonal)

Clone no. DF5

MONOSAN

Product name	Mouse anti-Beta1-Integrin, clone DF5 (Monoclonal)
Host	Mouse
Applications	IHC-P, IHC-fr, WB (1:1000), ELISA
Species reactivity	human
Conjugate	-
Immunogen	Wheat germ agglutinin binding glycoproteins of human fibroblasts
Isotype	IgG1
Clonality	Monoclonal
Clone number	DF5
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 1% BSA & 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Beta1-Integrin, clone DF5 (Monoclonal)

Clone no. DF5

MONOSAN

Additional info

The antibody reacts with the β_1 subunit of the integrin protein family and seems to be human specific. The antibody reacts with an extracellular epitope of the β_1 integrin molecule. Mab DF5 does react with paraffin sections

References

1. Koukoulis et al. Human Pathol 1997;28:84-90
2. Bernardini et al. 1997;83:673-678
3. Lohi et al. Histochem 1997;29:337-345
4. Virtanen et al. Am J Pathol 1995;150:1421-1431
5. Korhonen et al. Lab Invest 1990;62:616-625

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES