

Mouse anti-Beta1-Integrin, clone DF7 (Monoclonal)

Clone no. DF7

MONOSAN

Product name	Mouse anti-Beta1-Integrin, clone DF7 (Monoclonal)
Host	Mouse
Applications	IHC-P, IHC-fr, WB (1:1000), ELISA
Species reactivity	human
Conjugate	-
Immunogen	Wheat germ agglutinin binding glycoproteins of human fibroblasts
Isotype	IgG1
Clonality	Monoclonal
Clone number	DF7
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 1% BSA & 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Beta1-Integrin, clone DF7 (Monoclonal)

Clone no. DF7

MONOSAN

Additional info

The antibody reacts with the β_1 subunit of the integrin protein family. The antibody reacts with an extracellular epitope of the β_1 integrin molecule.

References

1. Ylänne et al. Int J Cancer 1989;43:1126-1136
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES