

Mouse anti-Chlamydia Trachomatis, clone 3 (Monoclonal)

Clone no. 3

MONOSAN

Product name	Mouse anti-Chlamydia Trachomatis, clone 3 (Monoclonal)
Host	Mouse
Applications	IHC-P, IF
Species reactivity	bacteria
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG3
Clonality	Monoclonal
Clone number	3
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Supernatant with 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Chlamydia Trachomatis, clone 3 (Monoclonal)

Clone no. 3

MONOSAN

Additional info

The antibody is directed to an immune dominant epitope (location unknown) of the Chlamydia LPS. This LPS is a genus specific antigen. As ascertained by the Dutch State Institute of National Health using a DOT-EIA the antibody reacts strongly with following prototype strains: - C. trachomatis: A B Ba C D E F G H I J K L L1 L2 L3 and mouse pneumonitis, -C. psittaci an ornithosis strain and a cat conjunctivitis strain, - C. pneumoniae: TW183

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES