

Mouse anti-alphaIIb-Integrin, clone CA3 (Monoclonal)

Clone no. CA3

MONOSAN

Product name	Mouse anti-alphaIIb-Integrin, clone CA3 (Monoclonal)
Host	Mouse
Applications	IHC-fr, WB (1:1000), ELISA
Species reactivity	human
Conjugate	-
Immunogen	Partially purified platelet glycoprotein IIb/IIIa complex
Isotype	IgG1
Clonality	Monoclonal
Clone number	CA3
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 1% BSA & 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti- α IIb-Integrin, clone CA3 (Monoclonal)

Clone no. CA3

MONOSAN

Additional info

The antibody reacts with the α II β -subunit of the β -integrins. Reacts specifically with platelets, megakaryocytes and leukemic cells able to megakaryocytic differentiation.

References

1. Tani et al. Exp Hematol 1996;24:158-168
2. Ylanne et al. Blood 1990;3:570-577
3. Ylanne et al. Blood 1989;72:1478-1486
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES