

Mouse anti-Cytokeratin 7,17,19, clone 4F5 (Monoclonal)

Clone no. 4F5

MONOSAN

Product name	Mouse anti-Cytokeratin 7,17,19, clone 4F5 (Monoclonal)
Host	Mouse
Applications	IHC-fr, WB (1:1000)
Species reactivity	human
Conjugate	-
Immunogen	Cytokeratins isolated from cultured pig kidney epithelial cells
Isotype	IgG1
Clonality	Monoclonal
Clone number	4F5
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 1% BSA & 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Cytokeratin 7,17,19, clone 4F5 (Monoclonal)

Clone no. 4F5

MONOSAN

Additional info

The antibody recognizes cytokeratins nos. 7, 17 and 19 typical for many simple epithelial cells but is found also in cells in squamous epithelia.

References

1. Rahman et al. Archs Dermatol Res 1993;285:17-31
2. Takumida te al. ORL 1991;53:10-14
3. Lehto et al. Cancer Rev 1986;6:1-20
4. Holthofer et al Lab Invest 1984;50:552
5. Virtanen et al. Ann N.Y. Acad.Sci 1984;455:635-648

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES