Product datasheet MON8062



Mouse anti-CD57, clone TB01 (Monoclonal)

Clone no. TB01 MONOSAN

Product name Mouse anti-CD57, clone TB01 (Monoclonal)

Host Mouse

Applications IHC-P,FC,IHC-fr

Species reactivity human

Conjugate Purified

Immunogen Human neuroblastoma cells.

lsotype lgM

Clonality Monoclonal

Clone number TB01

Size 0.1 mg

Concentration 1mg/ml

Format -

Storage buffer PBS with 0.09% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON8062



Mouse anti-CD57, clone TB01 (Monoclonal)

Clone no. TB01 MONOSAN

Additional info

Mouse anti Human CD57 antibody, clone TB01 recognizes CD57, also known as HNK-1, an oligosaccharide antigenic determinant present on a variety of polypeptides, lipids and chondroitin sulphate proteoglycans. Its function is poorly understood. CD57 is present on a subset of NK and T cells.

References 1. -

2 -

3. -

4. -5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES