

## Mouse anti-CD40, clone LOB7/6 (Monoclonal)

Clone no. LOB7/6

MONOSAN

---

Product name	Mouse anti-CD40, clone LOB7/6 (Monoclonal)
Host	Mouse
Applications	IHC-P,FC,IHC-fr,IP
Species reactivity	human,dog
Conjugate	Purified
Immunogen	CD40 Ig(Fc) fusion protein containing the EC region of human CD40 and Fc region of human IgG.
Isotype	IgG2a
Clonality	Monoclonal
Clone number	LOB7/6
Size	0.2 mg
Concentration	1.0 mg/ml
Format	-
Storage buffer	TBS
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD40, clone LOB7/6 (Monoclonal)

Clone no. LOB7/6

MONOSAN

**Additional info**

Mouse anti Human CD40 antibody, clone LOB7/6 recognizes the human CD40 cell surface antigen, a 48kDa glycoprotein expressed by B lymphocytes and weakly by some monocytes. CD40 is involved in the process of B cell selection in germinal centres and is vital in T cell-B cell interactions.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES