

## Mouse anti-CD21, clone 2G9 (Monoclonal)

Clone no. 2G9

MONOSAN

---

Product name	Mouse anti-CD21, clone 2G9 (Monoclonal)
Host	Mouse
Applications	IHC-P,IHC-fr
Species reactivity	human
Conjugate	S/N
Immunogen	Recombinant human CD21.
Isotype	IgG2a
Clonality	Monoclonal
Clone number	2G9
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD21, clone 2G9 (Monoclonal)

Clone no. 2G9

MONOSAN

**Additional info**

Mouse anti Human CD21 antibody, clone 2G9 recognizes the human CD21 cell surface antigen, a ~140 kDa glycoprotein. CD21 is expressed by mature B lymphocytes.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES