

Mouse anti-CD68, clone 514H12 (Monoclonal)

Clone no. 514H12

MONOSAN

Product name	Mouse anti-CD68, clone 514H12 (Monoclonal)
Host	Mouse
Applications	IHC-P,IHC-fr,IF
Species reactivity	human
Conjugate	S/N
Immunogen	Fusion protein corresponding to external domain of human CD68.
Isotype	IgG2a
Clonality	Monoclonal
Clone number	514H12
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD68, clone 514H12 (Monoclonal)

Clone no. 514H12

MONOSAN

Additional info

Mouse anti Human CD68 antibody, clone 514H12 recognizes the human CD68 cell surface antigen, a ~110 kDa glycoprotein primarily expressed by macrophages and monocytes.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES