Product datasheet MON8086



Mouse anti-Mcm2, clone CRCT2.1 (Monoclonal)

Clone no. CRCT2.1 MONOSAN

Product name Mouse anti-Mcm2, clone CRCT2.1 (Monoclonal)

Host Mouse

Applications WB,IHC-P

Species reactivity human

Conjugate Purified

Immunogen Bacterially expressed human Mcm2.

lsotype lgG1

Clonality Monoclonal

Clone number CRCT2.1

Size 0.1 mg

Concentration 1.0mg/ml

Format -

Storage buffer PBS with 0.09% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON8086



Mouse anti-Mcm2, clone CRCT2.1 (Monoclonal)

Clone no. CRCT2.1 MONOSAN

Additional info

Mouse anti Human MCM2 antibody, clone CRCT2.1 recognizes human DNA replication licensing factor MCM2, also known as Mcm2, Nuclear protein BM28 or mini chromosome maintenance protein-2. Mcm-2 is a 904 amino acid ~125kDa nuclear protein involved in the control of DNA replication. Mouse anti Human MCM2 antibody, clone CRCT2.1 has been used successfully for the detection of human MCM2 in ovarian adenocarcinoma by immunohistochemistry on formalin fixed, paraffin embedded tissues (Gakiopoulou et al. 2007).

References 1. -

2 -

3. -

4.

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES