

## Mouse anti-TGF beta, clone TB21 (Monoclonal)

Clone no. TB21

MONOSAN

---

Product name	Mouse anti-TGF beta, clone TB21 (Monoclonal)
Host	Mouse
Applications	WB,IHC-P,ELISA,FC,IHC-fr
Species reactivity	human,mink,sheep,rabbit,rat,mustelid,mouse
Conjugate	Purified
Immunogen	Human Transforming Growth Factor Beta 1 from platelets.
Isotype	IgG1
Clonality	Monoclonal
Clone number	TB21
Size	0.1 mg
Concentration	1.0 mg/ml
Format	-
Storage buffer	PBS
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Mouse anti-TGF beta, clone TB21 (Monoclonal)

Clone no. TB21

MONOSAN

**Additional info**

Mouse anti Human TGF beta antibody, clone TB21 recognizes both human platelet-derived and recombinant TGF-beta1 in enzyme-linked immunosorbent assay (ELISA). Mouse anti Human TGF beta antibody, clone TB21 demonstrates neutralising activity against TGF-beta1 in cell proliferation assays. Mouse anti Human TGF beta antibody, clone TB21 has been demonstrated to react with dimeric (&#126;25 kDa) or monomeric (&#126;12.5 kDa) molecules of natural TGF-beta1 under non-reducing and reducing conditions respectively.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES