

Mouse anti-Renal Cell Carcinoma, clone RC38 (Monoclonal)

Clone no. RC38

MONOSAN

Product name	Mouse anti-Renal Cell Carcinoma, clone RC38 (Monoclonal)
Host	Mouse
Applications	IHC-Fr (1:20)
Species reactivity	human
Conjugate	-
Immunogen	Renal cell carcinoma homogenate
Isotype	IgG1
Clonality	Monoclonal
Clone number	RC38
Size	1 ml
Concentration	15 ug/ml
Format	-
Storage buffer	Culture medium with 0.01M Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Renal Cell Carcinoma, clone RC38 (Monoclonal)

Clone no. RC38

MONOSAN

Additional info

The antibody stains over 95% of primary renal cell carcinomas and 60% of metastases of renal cell carcinoma and glomerular visceral epithelium and proximal tubules in normal kidney (cytoplasmic). It does not stain epithelial cells of non-renal malignancies.

References

1. Oosterwijk E et al. (1986). Am J Pathol 123, 301-309
2. Dijk J. van (1988) J Nucl Med 29, 1087-1082
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES