

Mouse anti-Carcinoembryonic Antigen, C.E.A., clone PARLAM4 (Monoclonal)

Clone no. PARLAM4

MONOSAN

Product name	Mouse anti-Carcinoembryonic Antigen, C.E.A., clone PARLAM4 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P (1:5-1:10), WB
Species reactivity	Human
Conjugate	-
Immunogen	Isolated from human colon carcinoma cells
Isotype	IgG1
Clonality	Monoclonal
Clone number	PARLAM4
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Carcinoembryonic Antigen, C.E.A., clone PARLAM4 (Monoclonal)

Clone no. PARLAM4

MONOSAN

Additional info

The antibody is directed against a single band of 180 kD of human carcinoembryonic antigen (CEA) and shows no cross-reactivity neither with bilary glycoprotein (BGP) nor with non-specific cross-reacting antigen (NCA).

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES