Product datasheet

MON9011P



MONOSAN

Mouse anti-MDR1 P-gp (ABCB1), clone JSB-1 (Monoclonal) Clone no. JSB-1

Product name	Mouse anti-MDR1 P-gp (ABCB1), clone JSB-1 (Monoclonal)
Host	Mouse
Applications	FC (1:10), IHC-fr, IHC-P (1:20), WB
Species reactivity	human, chinese hamster
Conjugate	-
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers
lsotype	lgG1
Clonality	Monoclonal
Clone number	JSB-1
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON9011P



Mouse anti-MDR1 P-gp (ABCB1), clone JSB-1 (Monoclonal)Clone no.JSB-1

MONOSAN

Additional info

The antibody reacts with a conserved cytoplasmic epitope of the plasma membrane-associated 170-180 kD glycoprotein, the expression of which is strongly correlated with the degree of multi-drug-resistance (MDR) derived MDR cell lines and human MDR cell lines, including cell lines derived from lung, ovaries and B cell lymphomas. Target species: Human, cross-reaction: Chinese hamster. No cross-reaction: mouse and rat.

Reference	c
VEIGI ELICE	2

1.

Scheper RJ et al. Int.J.Cancer 1988; 42: 389

- 2 Dalton WS et al. J Clin Oncol. 1989; 7: 415
- 3. Tiirikainen M et al. Ann. Haematol 1992; 65: 124
- 4. Grogan T et al. Blood 1993; 81: 490
- 5. Itsubo M et al. Cancer 1994; 73: 298

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com