

Mouse anti-Synaptophysin, clone SY38 (Monoclonal)

Clone no. SY38

MONOSAN

Product name	Mouse anti-Synaptophysin, clone SY38 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P (1:10), IF, WB
Species reactivity	Bovine, Human, Mouse, Rat
Conjugate	-
Immunogen	Synaptophysin from presynaptic vesicles, prepared from bovine brain
Isotype	IgG1
Clonality	Monoclonal
Clone number	SY38
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.5% BSA and 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Synaptophysin, clone SY38 (Monoclonal)

Clone no. SY38

MONOSAN

Additional info

The antibody is specific for Synaptophysin. The specificity was ascertained by immunoblotting and immunohistochemistry. The antibody predominantly reacts with a 38 kDa transmembrane glycoprotein from synaptic vesicles.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES