Product datasheet

MON9034



MONOSAN

Rat anti-MRP5 (ABCC5), clone M5II-54 (Monoclonal) Clone no. M5II-54

Product name	Rat anti-MRP5 (ABCC5), clone M5II-54 (Monoclonal)
Host	Rat
Applications	ICC (1:20-1:50), IHC-fr (1:20), FC, WB (1:50)
Species reactivity	human
Conjugate	-
Immunogen	MRP5 aa 722-910
lsotype	lgG2a
Clonality	Monoclonal
Clone number	M5II-54
Size	1 ml
Concentration	250 ug/ ml
Format	-
Storage buffer	supernatant with 0.7% BSA and 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON9034



Rat anti-MRP5 (ABCC5), clone M5II-54 (Monoclonal) Clone no. M5II-54

MONOSAN

Additional info

M5II-54 reacts with an internal epitope of MRP5, a 190-200 kD transmembrane protein that is closely related to the multidrug resistance protein MRP. M5II-54 was raised against a bacterial fusion protein of MRP5, containing amino acids 722-910 of the protein. M5I-1 does not cross-react with the human MDR1, MRP1, MRP2 or MRP3 gene products

References	1.	Kool et al. Cancer Res 1997; 57: 3537-3547
	2	Scheffer et al. Cancer Res 2000; 60: 5269-5277
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021