

Rat anti-MRP4 (ABCC4), clone M4I-10 (Monoclonal)

Clone no. M4I-10

MONOSAN

Product name	Rat anti-MRP4 (ABCC4), clone M4I-10 (Monoclonal)
Host	Rat
Applications	ICC, IHC-fr (1:20), WB (1:50), FC
Species reactivity	human, mouse
Conjugate	-
Immunogen	Maltose Binding Protein and aa 372-431 of MRP4
Isotype	IgG2a
Clonality	Monoclonal
Clone number	M4I-10
Size	1 ml
Concentration	100 ug/ ml
Format	Protein G purified
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-MRP4 (ABCC4), clone M4I-10 (Monoclonal)

Clone no. M4I-10

MONOSAN

Additional info

MRP4 transports cyclic nucleotides and anti-retroviral compounds. The M4I-10 Mab also reacts with the mouse orthologue of the transporter molecule (Mrp4).

References

1. Leggas M et al. Mol Cell Biol 2004; 24: 7612-7621
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES