Product datasheet MON9062-1



Mouse anti-CD120a, TNF-RI (p55/p60), clone MR1-2 (without preservatives) (Monoclonal)

Clone no. MR1-2 MONOSAN

Product name Mouse anti-CD120a, TNF-RI (p55/p60), clone MR1-2 (without preservatives)

(Monoclonal)

Host Mouse

Applications -

Species reactivity -

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

Clonality Monoclonal

Clone number MR1-2

Size 1 ml

Concentration n/a

Format -

Storage buffer PBS with 0.1% BSA and 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet MON9062-1



Mouse anti-CD120a, TNF-RI (p55/p60), clone MR1-2 (without preservatives) (Monoclonal)

Clone no. MR1-2 MONOSAN

Additional info

References 1.

2 -

3.

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES