

Mouse anti-CD120a, TNF-RI (p55/p60), clone MR1-2, biotinylated (Monoclonal)

Clone no. MR1-2

MONOSAN

Product name	Mouse anti-CD120a, TNF-RI (p55/p60), clone MR1-2, biotinylated (Monoclonal)
Host	Mouse
Applications	-
Species reactivity	-
Conjugate	Biotin
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Monoclonal
Clone number	MR1-2
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD120a, TNF-RI (p55/p60), clone MR1-2, biotinylated (Monoclonal)

Clone no. MR1-2

MONOSAN

Additional info

-

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES