

Mouse anti-Leukocyte Common Antigen, CD45 (Monoclonal)

Clone no. BRA-55

MONOSAN

Product name	Mouse anti-Leukocyte Common Antigen, CD45 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P (1:10-1:20), ELISA, IF, WB
Species reactivity	Human
Conjugate	-
Immunogen	Isolated from REH6 cell line
Isotype	IgG1
Clonality	Monoclonal
Clone number	BRA-55
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Leukocyte Common Antigen, CD45 (Monoclonal)

Clone no. BRA-55

MONOSAN

Additional info

The antibody reacts with 170 - 220 kD cell surface glycoproteins (CD45), the Leukocyte Common Antigen (LCA), expressed selectively on hematopoietic cells. It gives a pronounced staining on formalin-fixed, paraffin-embedded cells, and can be used to differentiate between lymphomas and carcinomas. Positive control: Tonsil.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES