

## Mouse ant-Collagenase G, clone 1D2 (Monoclonal)

Clone no. 1D2

MONOSAN

---

Product name	Mouse ant-Collagenase G, clone 1D2 (Monoclonal)
Host	Mouse
Applications	ELISA(1:500-1:5000), WB(1:2000-1:5000)
Species reactivity	human
Conjugate	-
Immunogen	Collagenase from Clostridium histolyticum
Isotype	IgG2b
Clonality	Monoclonal
Clone number	1D2
Size	100 ug
Concentration	0.75 mg/ml
Format	-
Storage buffer	PBS, pH 7.2
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse ant-Collagenase G, clone 1D2 (Monoclonal)

Clone no. 1D2

MONOSAN

**Additional info**

This Collagenase is typically used with a neutral protease for tissue disruption.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES