

Mouse anti-EBFP, clone 3A6 (Monoclonal)

Clone no. 3A6

MONOSAN

Product name	Mouse anti-EBFP, clone 3A6 (Monoclonal)
Host	Mouse
Applications	ELISA, WB, IP
Species reactivity	human
Conjugate	-
Immunogen	Recombinant EBFP (NCBI accession number AX_766758 REGION: 1-717)
Isotype	IgG1
Clonality	Monoclonal
Clone number	3A6
Size	100 ug
Concentration	n/a
Format	-
Storage buffer	lyophilized in a 10 mM ammonium bicarbonate buffer. Each vial contains 2 mg BSA
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-EBFP, clone 3A6 (Monoclonal)

Clone no. 3A6

MONOSAN

Additional info

Fluorescent proteins, like EBFP, can be used as protein "tags" to study the subcellular localization of proteins and/or their translocation upon stimulation and/or as markers for transfections in transient and stable expression systems. In immunoblot it cross reacts with GFP and YFP.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES