

Mouse anti-Prostate Specific Antigen (PSA), clone ER-PR8 (Monoclonal)

Clone no. ER-PR8

MONOSAN

Product name	Mouse anti-Prostate Specific Antigen (PSA), clone ER-PR8 (Monoclonal)
Host	Mouse
Applications	IHC-fr, IHC-P (1:10), IF, WB
Species reactivity	Human
Conjugate	-
Immunogen	Isolated from human prostatic carcinoma
Isotype	IgG1
Clonality	Monoclonal
Clone number	ER-PR8
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Prostate Specific Antigen (PSA), clone ER-PR8 (Monoclonal)

Clone no. ER-PR8

MONOSAN

Additional info

The presence of prostate specific antigen is highly specific for demonstration of the prostatic origin of a malignancy. Mab ER-PR8 has been shown to be a valuable reagent for the identification of primary and metastatic prostatic carcinoma. The antibody reacts with a 34kD protein in benign and malignant prostatic epithelium.

Positive control: Prostatic benign hyperplasia.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES