Product datasheet MON9128



Rabbit anti-Human NKX3.1, clone EP356

Clone no. EP356 MONOSAN

Product name Rabbit anti-Human NKX3.1, clone EP356

Host Rabbit

Applications IHC-P (1:50-1:200)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

lsotype lgG

Clonality Monoclonal

Clone number EP356

Size 1 ml

Concentration n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES



Rabbit anti-Human NKX3.1, clone EP356

Clone no. EP356 MONOSAN

Additional info

NKX3.1 is a prostate specific androgen-regulated homeobox gene located on chromosome 8p. It is difficult to distinguish between high grade prostate adenocarcinoma and high grade infiltrating urothelial carcinoma using hematoxylin and eosin stained specimens. Current prostate adenocarcinoma markers such as prostate specific antigen (PSA) and prostate specific acid phosphatase (PSAP) are very useful in determining prostate origin of prostate cancer in other sites, but have lower sensitivity when identifying poorly differentiated compared to well differentiated cases. NKX3.1 is a sensitive and specific tissue marker of prostate adenocarcinoma and can be used to help distinguish it from urothelial carcinomas. Currently, thrombomodulin and uroplakin are used to identify tumors of urothelial origin; however, their sensitivities are suboptimal.

References 1. Gurel B, et al. Am | Surg Pathol. 2010; 34:1097-105

2 Chuang AY, et al. Am J Surg Pathol. 2007; 31:1246-55

3. -

4. -

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES