

Rat anti-Mouse MRP6, clone M6II-24 (Monoclonal)

Clone no. M6II-24

MONOSAN

---

Product name	Rat anti-Mouse MRP6, clone M6II-24 (Monoclonal)
Host	Rat
Applications	ICC, IHC-fr (1:20-1:50), WB
Species reactivity	mouse
Conjugate	-
Immunogen	mouse MRP6 aa 843-999
Isotype	IgG2a
Clonality	Monoclonal
Clone number	M6II-24
Size	1 ml
Concentration	50 ug/ ml
Format	-
Storage buffer	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse MRP6, clone M6II-24 (Monoclonal)

Clone no. M6II-24

MONOSAN

**Additional info**

M6II-24 reacts with mouse Mrp6, a 165 kD transmembrane protein that is related to the multidrug resistance related protein Mrp1. Mutations in the human MRP6 gene are responsible for the connective tissue disorder Pseudoxanthoma elasticum (PXE). M6II-24 was raised against a bacterial fusion protein of mouse Mrp6, containing amino acids 843-999.

**References**

1. Bergen AA et al. Nat Genet 2000; 25: 228-231
2. Scheffer GL et al. Lab Invest 2002; 82: 515-518
3. Gorgels TG et al. Hum Mol Genet 2005; 14: 1763-73
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**