### Product datasheet

MON9825



MONOSAN

Mouse anti-Blood Group H Type 5,6 (CD173), clone HE-182 (Monoclonal)Clone no.HE-182

Product name	Mouse anti-Blood Group H Type 5,6 (CD173), clone HE-182 (Monoclonal)
Host	Mouse
Applications	ELISA, FC, IHC-fr, IHC-P, IF
Species reactivity	human
Conjugate	-
Immunogen	human erythrocytes
lsotype	IgM-K
Clonality	Monoclonal
Clone number	HE-182
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

# FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

### Product datasheet

MON9825



Mouse anti-Blood Group H Type 5,6 (CD173), clone HE-182 (Monoclonal)Clone no.HE-182

MONOSAN

#### Additional info

HE-182 reacts with human blood group H type 5 and type 6 as confirmed with Chembiomed synthetic oligosaccharide-BSA conjugate in ELISA. It is expressed on endothelial cells, epithelial cells and granulocytes. Increased expression of this antigen has been observed on some tumor tissues such as gastric carcinomas, urothelial carcinomas, and colon carcinomas

References	1.	Rouher Ph et al. Blood transfusion and immuohaematology, Arnette France 30
	2	-
	3.	-
	4.	-
	5.	-

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021